

#### AREA ARACE ORACE O

PUBLIC ACCESS BLEEDING CONTROL KITS ANSI WORKPLACE FIRST AID TRAUMA KITS

TRAUMA & FIRST AID KIT



The first person on scene IS the first responder...

> Are you Prepared?

Kit shown: Class B Trauma & First Aid Kit-Hardcase see page 4

ATT THE AVE

## The Difference between Life and Death

**MINUTES** 

3to5

A victim who is bleeding can quickly die BEFORE professional emergency responders can even arrive on the scene

- Serious bleeding accounts for 35% of prehospital deaths from traumatic injuries
- Life-threatening bleeding warrants immediate intervention
- Basic bleeding control provided by bystanders can often be the difference between life and death...

## **FACT:** from 2000-2020, **45**% of ALL Active Shooter Incidents occurred in Places of Commerce





Uncontrolled BLEEDING can occur in as little as 3 to 5 Minutes

Trauma is a national health problem and after a traumatic injury, serious bleeding is responsible for over 35% of prehospital deaths.\*\* Causes of trauma can range from natural disasters, falls, traffic, industrial, and workplace accidents to intentional acts of violence such as shootings, stabbings and bombings. In fact, according to data collected from 2000 to 2020. the U.S. Department of Justice found that the workplace is the most likely location of an active shooter situation.\*

Because these traumatic injuries often result in uncontrolled bleeding, immediate intervention – even before emergency medical services or other first responders can arrive – is critical to avoid potential death.

**Death from** 

#### North American Rescue

is leading the way to equip bystanders and initial first responders with the essential medical gear needed to take immediate action in life threatening emergencies. These lifesaving modular aid components are designed for individuals or groups in the workplace, public buildings, and populated public areas to administer immediate bleeding control at the point of wounding.

\*Active Shooter Incidents 20-Year Review, 2000-2020, Federal Bureau of Investigation, U.S. Department of Justice, Washington, D.C.

\*\*Chambers JA, Seastedt K, Krell R, Caterson E, Levy M, Turner N. "Stop the Bleed": A U.S. Military Installation's Model for Implementation of a Rapid Hemorrhage Control Program. Mil Med. 2019 Mar 01;184(3-4):67-71. [PubMed] [Reference list] https://academic.oup.com/milmed/article/184/3-4/67/5062295

## Workplace First Aid Kits

# Class A TRAUMA & FIRST AID KIT

- Designed as a portable First Aid Kit for the workplace exceeding ANSI/ISEA Z308.1-2021 Standard for Class A kits
- Able to treat the most common workplace injuries with the added benefit of being able to control major bleeding with the C-A-T Tourniquet the safest and most effective Tourniquet on the market!
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

#### **TRAUMA & FIRST AID KIT - CLASS A**

CASE	DRESSING	ITEM #
Soft Case	Compressed Gauze	900060823
Soft Case	Bleeding Control Dressing	900060824
Hard Case	Compressed Gauze	900060825
Hard Case	Bleeding Control Dressing	900060826

#### SPECIFICATIONS

**Soft Case:** L 14.75 in. x H 8 in. x D 4 in. Weight: 2.94 lb

TRAUM

COLDSTAR

USA

TRAUMA & First aid kit

Hard Case: L 13 in. x H 9.5 in. x D 5 in. Weight: 6 lb

#### TRAUMA & FIRST AID KIT CLASS A:



## Workplace First Aid Kits

## Class B TRAUMA & **FIRST AID KIT**

- Designed as a portable First Aid Kit for the workplace exceeding ANSI/ISEA Z308.1-2021 Standard for Class B kits
- Same supplies as the Class A Kit but with larger quantities to accommodate larger workplaces
- Packaged in a convenient zippered nylon bag or rugged hard case that is easily transported to the patient's side

0

#### **TRAUMA & FIRST AID KIT - CLASS B**

SPECIFICATIONS Soft Case:	CASE	DRESSING	ITEM #	
_ 14.75 in. x H 8 in. x D 7.5 in.	Soft Case	Compressed Gauze	900060827	
Weight: 4.87 lb	Soft Case	Bleeding Control Dressing	900060828	
<b>Hard Case:</b> L 13.5 in. x H 11 in. x D 6 in.	Hard Case	Compressed Gauze	900060829	
Weight: 9 lb	Hard Case	Bleeding Control Dressing	900060830	

#### IN TRAUMA & FIRST AID KIT CLASS B:



### Kits & Stations | Public Access Bleeding Control

**SPECIFICATIONS** 

Nylon Bag: H 7 in. x W 5 in. x D 2.5 in.

Vacuum Sealed: H 7 in. x W 5 in. x D 2.75 in.

## Public Access INDIVIDUAL **BLEEDING CONTROL KITS**

#### Nylon Bag Option:

- · Easy to restock in cases of use
- Convenient to carry and easy to open in stressful situations
- Opens flat for ease of use and immediate identification of enclosed devices
- Easily recognizable red nylon bag with attached PVC Bleeding Control ID Patch

#### Vacuum Sealed Option:

- Cost-effective alternative to Nylon
- Compact and rugged Vacuum-Sealed pouch
- Large Identifying Card Vacuum-Sealed inside the pouch for ease of identification

#### Every INDIVIDUAL BLEEDING CONTROL KIT contains:



#### 8-PACK CLEAR POLYCARBONATE WALL STATIONS

Gauze<sup>®</sup>

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	8	Vacuum	21.8 lb	900060819
Advanced BCD	8	Vacuum	23 lb	900060834
Basic	8	Nylon	23.5 lb	900060835
Advanced BCD	8	Nylon	27.4 lb	900060836



#### INDIVIDUAL BLEEDING CONTROL KITS

LEVEL	BAG	WEIGHT	ITEM #
Basic	Red Nylon	15 oz	900060817
Intermediate	Red Nylon	16.5 oz	900060820
Advanced BCD	Red Nylon	17.4 oz	900060831
Basic	Vacuum Sealed	11.7 oz	900060816
Advanced BCD	Vacuum Sealed	14 oz	900060822

## Clear Polycarbonate WALL STATIONS

- Includes Individual Bleeding Control Kits packaged in red nylon bags or Vacuum-Sealed pouches.
- The 8 Pack Stations also include an easy to-open nylon carrier and one (1) Responder QuikLitter™ for moving injured casualties.



#### **SPECIFICATIONS**

8 Pack Stations: H 18 in. x W 15.5 in. X D 10 in. Empty Case Wt: 15 lb (includes rack) Individual Stations: H 10 in. x W 10 in. X D 5.25 in. Empty Case Wt: 2.39 lb

#### **Individual Kit Station**



#### IND. CLEAR POLYCARBONATE WALL STATIONS

LEVEL	#KITS	BAG	WGT	ITEM #
Basic	1	Vacuum	3.3 lb	900060837
Advanced BCD	1	Vacuum	3.48 lb	900060838
Basic	1	Nylon	3.3 lb	900060839
Advanced BCD	1	Nylon	3.7 lb	900060840

## Individual Components

## C-A-T<sup>®</sup> Combat Application Tourniquet Gen 7

- Single route buckle decreases application time to help decrease blood loss, allows for slack removal for proper TQ application
- Applied same way to both upper and lower extremities
- True one-handed tourniquet

.....

• Scientifically proven 100% effective in occluding blood flow in both upper and lower extremities by the USAISR and recommended by CoTCCC

#### QuikClot BLEEDING **CONTROL DRESSING**

- Advanced hemostatic agent combined with the simplicity of gauze available without medical device authorization (MDA)
- Utilizes the same inorganic component as QuikClot Combat Gauze<sup>®</sup> used by the US Military to control severe bleeding

#### QuikClot Bleeding Contro! Dressing

For Temporary External Vie to Control Bleeding um Packed and Z-Folded

#### SPECIFICATIONS L 4 yd x W 3 in. Pkg: H 5 in. x W 5.5 in. x

## D 0.625 in.

Weight: 0.8 oz

#### **QUIKCLOT BCD**

ITEM #

900060821

#### Emergency Trauma Dressing **RESPONDER ETD**

- Suitable for applying immediate direct pressure
- Equipped with a non-adherent pad
- Easy-to-use hook & loop securing device with no confusing pressure bars or hooks
- Simple application with the Quick grip roll control system which prevents the dressing from unrolling if dropped.
- Rugged vacuum-sealed, low-cube packaging with easy-open tear notches

#### **RESPONDER EMERGENCY TRAUMA DRESSING**

BANDAGE	PAD	WGT	PACKAGE	ITEM #
W6 in. x L60 in.	6x6	2.9 oz	L6 in. x W2.5 in. x D2.1 in.	900060832

#### NAR WOUND PACKING GAUZE

- Ultra compressed, soft, white 1-ply gauze includes embedded X-Ray detectable line
- Z-fold packaging for controlled application as a complete unit or in a continuous linear feed
- Used as a packing gauze for wounds, as a backing gauze for hemostatic agents, with a compression bandage to control hemorrhage, or for minor wound bandaging

#### **SPECIFICATIONS**

Gauze: W 4.5 in. x L 4.1 yd Packaged: L 3 in. x W 2 in. x D 1 in. Weight: 1.2 oz

#### WOUND **PACKING GAUZE**

ITEM #

900060818

#### SURVIVAL BLANKET Aluminized, reflective PET to help retain up to 90% of your body heat

**C-A-T TOURNIQUET** 

COLOR

Orange

Black

Blue

NAR

STYLE

Tactical

Rescue

Trainer

# IIIIIIII

**SPECIFICATIONS** 

**ITEM #** 

900060812

900060813

900060815

Device: L 37.5 in.

Weight: 2.7 oz

H 6.5 in. x W 2.4 in. x D 1.5 in.

Packaged:

#### **SPECIFICATIONS**

Unfolded: L 84 in. x W 52 in. Weight: 1.7 oz

**ITEM #** NSN # 900060833 4240-01-697-8869

**HYFIN® VENT COMPACT Twin Pack** 

- Two Chest Seals packaged approximately 25% smaller than our standard HyFin®
- 3-channel pressure relief vents
- Two Chest Seals designed for the treatment of entry/exit or multiple penetrating wounds
- Large red tab designed for single step, peel and apply application and for burping the wound if necessary

#### **SPECIFICATIONS**

Allin A

Device (excludes tab) 4.75 in. x 4.75 in. Wqt: 1.55 oz Packaged: H 3.88 in x W 6.6 in x D 0.13 in. U.S. Patent: 7,504,549

#### **HYFIN VENT COMPACT ITEM #** 900060814

www.lawsonproducts.com • 800-890-8198







www.lawsonproducts.com Government Customer Care 800-890-8198 government.sales@lawsonproducts.com

ASS

HIT



CPR, AED, and

PUBLIC ACCESS BLEEDING CONTROL KITS | ANSI WORKPLACE FIRST AID TRAUMA KITS

REV021623